



ANESTHESIA MONITORING QUICK REFERENCE

PARAMETER	RANGE
Heart rate (anesthetized)	Dogs: 60–140 bpm* Cats: 100–200 bpm*
Respiratory rate (anesthetized)	8–20 bpm*
Temperature (anesthetized)	99.0–101.5° F 37.2–38.6° C
Systolic arterial pressure (SAP)	90–140 mmHg*
Mean arterial pressure (MAP)	65–90 mmHg*
Diastolic arterial pressure (DAP)	55–90 mmHg*
SpO ₂	95–100%
ETCO ₂	35–45 mmHg
FiCO ₂	0–5 mmHg
Isoflurane MAC	Dogs: 1.25% Cats: 1.6%
Sevoflurane MAC	Dogs: 2.4% Cats: 2.6%

*Values may differ depending on the specific medications administered.

Interested in anesthesia and patient safety training for your hospital? Thomson Veterinary Education develops structured, practical, and cost-effective training programs and institutional safety tools to meet your hospital's needs. Visit [ThomsonVeterinaryEducation.com](https://www.ThomsonVeterinaryEducation.com) for more details.